



Get Breast Milk Home from Anywhere in the U.S.

portal.milkstork.com/exelixis




How It Works




PUMP

Pump according to your regular nursing schedule and refrigerate or freeze your breast milk.



PACK

Pack your milk into the cooler provided. If shipping, your cooler will be pre-labeled with a FedEx Priority Overnight shipping label.



SHIP

Ship your milk from an authorized FedEx shipping location or schedule a FedEx pick-up with your hotel or destination office.




TOTE/CHECK

Tote your cooler or check it as luggage using the Milk Stork travel bags provided.

Let Milk Stork Do the Heavy Lifting




Simple Order Process & Travel Logistics



Pre-Labeled Ready-To-Go Refrigerated & Frozen Kits



Expedited Shipping on All Deliveries



Proactive Shipment Monitoring

Need more info?
Visit support.milkstork.com